

# Longest™

## Compression Therapy System DVT Sleeve



Longest™ Compression Therapy System DVT Sleeve (single patient use) as an essential part of Longest™ Compression Therapy System is compatible with both Longest™ Compression Therapy Device LGT-2200DVT and LGT-2201DVT.

It is made of premium fabric and TPU with featuring lightweight, breathable, comfortable which is able to help prevent a build-up of heat and moisture contributing to increase patient's compliance during therapy.

With an extensive range of foot, calf, thigh, foot and calf, thigh and calf sleeves, Longest™ DVT Sleeve address a wide range of patient types, and caregiver needs and preferences.

## ⊕ Feature

- ⊗ Made of premium fabric with more lightweight, breathable and comfortable
- ⊗ Compatible with both Longest™ LGT-2200DVT and LGT-2201DVT
- ⊗ With a wide range of types for different application
- ⊗ One-Three (1-3) chambered sleeve with intermittent and sequential (gradient) inflation

## DVT Sleeve Size Details

| Order No. | Sleeve               | Size | Proximal Circumference |
|-----------|----------------------|------|------------------------|
| 220092    | Thigh Length Sleeve  | S    | <55.9cm                |
| 220092M   |                      | M    | <71.1cm                |
| 220092L   |                      | L    | <91.4cm                |
| 220093M   | Calf and Foot Sleeve | M    | <53.3cm                |
| 220094    | Calf Sleeve          | S    | <40.6cm                |
| 220094M   |                      | M    | <48.3cm                |
| 220094L   |                      | L    | <55.9cm                |
| 220095M   | Foot Sleeve          | M    | <40.6cm                |



**Longest™**

**Guangzhou Longest Medical Technology Co., Ltd.**

Address: #96 Chuangqiang Road, Ningxi Street, Zengcheng District, Guangzhou, China 511399

Tel: 86-20-66353999 Fax: 86-20-66353920

E-mail: export@longest.cn Website: www.longestmedical.com

