

Certification



ISO13485 / MDSAP



Health Canada



CE

Specification

Output Channel	Ⓐ 12 Channels consist by two groups; Ⓑ Each channel can be set and adjusted independently.
Treatment Mode	Nine Modes Standard, Toning, Fast Slim, Lymph Drainage, Relaxation, Abdominals, Thighs&Hips, Facial, Bust
Output Intensity	Adjustable, 0-120mA
Treatment Time	1-99min, stepping 1min
Speed (Training frequency)	Slow, Medium and Fast
Safety Protection	Over-current Protection Over-voltage Protection Short Circuit Protection

Longest

Guangzhou Longest Science & Technology Co., Ltd.
Address: #96 Chuangqiang Road, Ningxi Street, Zengcheng District,
Guangzhou, China 511399
Tel: 86-20-66353999 Fax: 86-20-66353920
E-mail: export@longest.cn Website: www.gzlongest.com



LGT-2320BE

Electrotherapy Device for Medical Aesthetics

Health Life Longest Care

Longest



- ⌄ Lose weight
- ⌄ Tone your body
- ⌄ Get rid of cellulite

Electronic Muscle Stimulation for Beauty


(Model Number:LGT-2320BE) applies low-frequency currents (modulation or non-modulation) to stimulate the motor nerves or muscles to make them contract. After the muscle contractions, the muscle fiber thickens, the weight of the muscle volume increases, the capillaries in the muscle becomes rich, which strengthens the entire muscle to achieve the effect of losing weight and shaping.



Features

- Multiple Treatment Areas
- No weights during your workout
- Non-Invasive, No Downtime
- A 45-minutes session can burn up to 1000 calories
- Intensity automatic locking for your safety

Application



- | | |
|-------------------|---|
| Back of the Chest | Relaxation Mode
Build up Muscle, Strengthening and Relaxing |
| Upper Arms | Standard Mode
Build up Muscle, Strengthening and Toning |
| Abdomen | Abdominals Mode
Cellulite Reduction, Strengthening and Toning |
| Buttocks | Thigh&Hip Mode
Fat Breakdown and Toning |
| Thigh | Thigh&Hip Mode
Break down fat, Strengthening and Lymphatic Drainage |
| Lower Leg | Toning Mode
Toning and Cellulite Reduction |

Accessory



Silicone Rubber Electrode*12



Electrode Lead Hose
(Φ95mm) *12



Straps with 3 Sizes*12
(4 pcs per size)



Safety Switch*1

