




Certification




Promoting member of ISMST




CE (MDD)




ISO13485 / ISO9001




Health Canada




Therapeutic Goods Administration



FDA



NRTL



INMETRO

Technical Specification

Pressure Range:	1.0~4.0 Bar (Step by 0.1 Bar)
Frequency:	1~17Hz
Pre-set Protocol:	Yes
Circumference:	Yes
Dimension:	39.5 (L) x 39 (W) x 19.4 (H) cm
Standard Weight:	11.3kg
Langages:	English, Italian, German, Spanish, French

Longest

Guangzhou Longest Science & Technology Co., Ltd.  
Address: #96 Chuangqiang Road, Ningxi Street, Zengcheng District,  
Guangzhou, China 511399  
Tel: 86-20-66353999    Fax: 86-20-66353920  
E-mail: export@longest.cn    Website: www.gzlongest.com



PowerWave  
LGT-2500S

The ideal body shaping treatment  
Acoustic waves provide an ideal body shaping treatment option

Health Life    Longest Care

Longest





# PowerWave LGT-2500S

## Acoustic Wave Therapy Device

Acoustic wave treatment uses acoustic sound waves together with the body's own ability to repair tissue by stimulating circulation in all layers of the skin, including blood vessels, connective tissue and muscles.

AWT tightens and forms new connective tissue, reducing cellulite and improving skin appearance.

## ⊕ Applications

- ▶ Anti-cellulite treatment
- ▶ Body reshaping
- ▶ Skin elasticity improvement
- ▶ Treatment of skin irregularity after liposuction
- ▶ Connective tissue tightening
- ▶ Scar and wrinkle smoothing
- ▶ Treatment of stretch marks
- ▶ Lymphatic system stimulation



## ⊕ Benefits

- ▶ The acoustic wave is able to massage the tissues in depth with considerable beneficial effects.
- ▶ Stimulation of blood circulation and lymphatic
- ▶ Raising the heat tissue and activation of the elimination of fat cells
- ▶ Activation of metabolism
- ▶ Stimulating the production of collagen, improving elasticity and texture of the skin



## ⊕ Effects

- ▶ Body contouring
- ▶ Cellulite treatment
- ▶ Firming tissue
- ▶ Smoothing of scars
- ▶ Improvement in skin elasticity
- ▶ Treatment of pregnancy stretch marks
- ▶ Circumference reduction

